

Carbon / Water Analyzer CW-800M General Information

ELTRA's CW-800M analyzer is designed for fractional analysis of carbon and water in one single operation. It is based on the proven technology of the CW-800 but has a modified furnace which allows fast adjustment of the temperature in the furnace. Different fractions of carbon and water of the sample are analyzed by applying different furnace temperatures. ELTRA's CW-800M analyzer is designed for the precise, simultaneous determination of carbon (released as CO2) and water in lime, gypsum and cement from trace level up to 100 % (depending on sample weight). Other sample materials include ores, soil, minerals, slags and waste.

Up to five temperature steps ("ramps") can be programmed for each application. The temperature levels and their durations are selectable. The maximum temperature is $1000\,^{\circ}$ C.

Depending on the application step, either O2 can be used as furnace atmosphere (oxidation of the sample) or inert gas like N2 or argon.

The detection system of ELTRA's CW-800M is very sensitive, reliable and guarantees a long lifetime. It can be customized according to the user's requirements. Two infrared cells can be combined independently and allow highly precise measurement of the released CO2 and H2O.



cement, gypsum, limestone, minerals, ores, slag, soil, waste

Product Advantages

- simultaneous carbon dioxide and water determination with minimal sample preparation
- Analysis of TOC (Total Organic Carbon) and TIC (Total Inorganic Carbon) without adding acids
- rapid, precise, accurate and reliable element determination
- · wide range of materials can be analyzed
- resistance furnace temperature can be set up to 1000 °C in steps of 1 °C
- up to 5 programmable steps with different temperatures ("ramps") can be defined
- · customized infrared cells provide wide, dynamic measuring range
- due to gold IR path, increased cell live time for analysis of halogen or acid containing samples
- powerful software (multilingual, customized display, export of results)
- single and multipoint calibration
- no halogen trap required
- · electronic gas flow control
- low maintenance
- robust design allows usage in production control and laboratory

Features





Carbon / Water Analyzer CW-800M

Measured elements carbon dioxide, water

Furnace alignment horizontal Sample carrier quartz boats

Field of application agriculture, biology, chemistry /

plastics, construction materials, environment / recycling, geology /

mining, others

Furnace resistance furnace with quartz tube,

adjustable up to 1000 °C

Catalyst furnace +

Process of measurement temperature and carrier gas can be

changed during measurement according to a user-defined program

Detection method solid state infrared absorption

Number of IR cells 1 - 2 Material of IR path gold

Typical analysis time 5 - 20 min (depending on program)

Chemicals required copper oxide, magnesium

perchlorate, sodium hydroxide

Gas required nitrogen 99.995 % pure (2 - 4 bar /

30 - 60 psi) oxygen 99.995 % pure

(2 - 4 bar / 30 - 60 psi)

Power requirements 230 V, 50/60 Hz, max. 10 A, 2300 W

Dimensions (W x H x D) 55 x 80 x 60 cm

Weight ~ 65 kg

Required equipment balance (resolution 0.0001g),

monitor, PC

Optional accessories TIC module, voltage stabilizer 5 KVA

Function Principle

Operation of the CW-800M is simple and convenient. After weighing the sample in a quartz boat, it is placed on the loading mechanism of the furnace. In the following, the analysis can be started and the boat is introduced into the furnace by the user. Depending on a user defined program, different temperatures and carrier gases are applied to the sample. While processing this program, the released CO2 and water is determined by the infrared cells. The received "chromatogram" of the sample subsequently shows different fractions of carbon and water of the sample. All data processing, control of the combustion process and calculating of the result is done by an external PC. The duration of the measurement depends on the length of the user defined steps. A common analysis takes about up to 20 minutes.

incl. order data

ELTRA CW-800M "Multiphase"



Carbon / Water Analyzer CW-800M

(Please order PC, monitor, balance and consumables (starter-kit, anhydrone, sodium hydroxide, copper oxide) separately)

Measuring ranges at 200 mg sample weight

88100-4042 CW-800M H2O 0 - 20 % H2O

88100-4041 CW-800M 1x CO2 0 - 70 % CO2 | + H2O 0 - 20 %

H2O

88100-4055 CW-800M 1x CO2 0 - 1 % CO2 | + 1x CO2 1 - 70 %

CO2

PC, Monitor, Balance

71015 Computer with dual core processor, 300 GB HDD, 4

GB RAM, Windows operating system, DVD-ROM,

keyboard, mouse

71016 Monitor, TFT

88600-0002 Balance (resolution 0.0001 g)

71002 Printer

Accessories

38001 TIC-Module
72070 Oxygen regulator
72080 Nitrogen regulator
71090 Voltage Stabilizer 5 KVA

Consumables

Required consumables

88500-0014 Starter-kit for 1,000 analyses (1,000 disposable

porcelain boats, 3 quartz boats, 50 g glass wool, 50 g

quartz wool)

90200 Anhydrone (magnesium perchlorate), 454 g

90210 Sodium hydroxide, 500 g 90290 Copper oxide, 100 g

Optional consumables

36120 Quartz boat 75x7.4x1.5 mm

90160 Disposable porcelain boats 86x13x10 mm, 1,000

pieces

90332 Glass wool, 50 g 90330 Quartz wool, 50 g

90810 Calcium carbonate, 100 g 90700-1040 Calcium oxalate, 50 g

92610 Tube of high vacuum grease



Carbon / Water Analyzer CW-800M

Spare and Wear Parts

77503 Heating element 48750 Combustion tube