



PACIFIC
BIOSCIENCES®



Technical Overview: Multiplexed Library Preparation for Full-Viral Genome Sequencing Using HiFiViral SARS-CoV-2 Kit

Sequel II and IIe Systems ICS v10.1 / SMRT Link v10.2

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2021 by Pacific Biosciences of California, Inc. All rights reserved.

PN 102-205-300 Version 01 (November 2021)

Multiplexed Library Preparation for Full-Viral Genome Sequencing Using HiFiViral SARS-CoV-2 Kit


1. HiFiViral SARS-CoV-2 Kit Workflow Overview
 2. Multiplexed Library Preparation Using Molecular Inversion Probe-Based Enrichment with the HiFiViral SARS-CoV-2 Kit
 3. Multiplexed SARS-CoV-2 Library Sequencing Workflow Recommendations
 4. Multiplexed SARS-CoV-2 Data Analysis Recommendations
 5. Multiplexed SARS-CoV-2 Library Example Performance Data
 6. Technical Documentation & Applications Support Resources
- APPENDIX 1:* RNA Isolation Kit Options for Full-Viral Genome Sequencing of SARS-CoV-2
- APPENDIX 2:* Guidance on Workflow Automation For Multiplexed Library SARS-CoV-2 Library Preparation

SARS-CoV-2 FULL-VIRAL GENOME SEQUENCING: HOW TO GET STARTED




[PacBio COVID-19 Sequencing Tools and Resources Website](#)

Summary overview of application-specific sample preparation and data analysis workflow recommendations




Procedure & Checklist: HiFiViral SARS-CoV-2 Workflow ([102-188-800](#))*

Technical documentation containing sample library construction and sequencing preparation protocol details



Technical Overview: Multiplexed Library Preparation for Full-Viral Genome Sequencing Using HiFiViral SARS-CoV-2 Kit ([102-205-300](#))

Technical Overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.



* For Research Use Only. Not for use in diagnostic procedures.

HiFiViral SARS-COV-2 KIT USES MOLECULAR INVERSION PROBES FOR EFFICIENT ENRICHMENT OF VIRAL RNA SEQUENCES FOR ANALYSIS



Robust performance



Easier workflow



Capture all variants



Flexible batch size



Cost-effective

Better Performance with Molecular Inversion Probes (MIPs)

- Differentiated enrichment technology
- Robust genome coverage across a range of Ct-values
- Probe design resilient to novel variants
- Capture mutations of all types
- Detect multiple strains in one sample


Easier Workflow and Faster Turnaround Times

- Easier workflow compared to targeted PCR amplicons
- All ready-to-use reagents in one kit
- Color change indicator confirms correct reagent was added
- Addition-only workflow can be automated
- Automated sequencing and analysis runs overnight


Flexible Scaling

- 384 reactions per kit
- Scalable batching: 24 – 384 samples per run

Quickly and efficiently scale genomic surveillance by sequencing with **an accurate and robust kit solution to capture all variants**



END-TO-END PACBIO PROTOCOL FOR FULL-VIRAL GENOME SEQUENCING USING HiFiViral SARS-CoV-2 KIT



HiFiViral SARS-CoV-2 Kit ([102-132-000](#))

For Targeted Enrichment and Barcoding of SARS-CoV-2 PCR-Positive Samples*



- ✓ Robust Performance
- ✓ Easier Workflow
- ✓ Capture All Variants
- ✓ Flexible Batch Size
- ✓ Cost Effect


PROCEDURE & CHECKLIST

PacBio HiFiViral High Throughput Multiplexing for Full-Viral Genome Sequencing of SARS-CoV-2

November 2021

PacBio Documentation ([102-188-800](#))

- Full workflow can be completed from **sample to answer in as short as ~28 – 42 h** (1 – 2.5 h hands-on time)
- Multiplex 24 – 384 samples per SMRT Cell 8M and load up to 8 SMRT Cells per Sequel IIe System to **run up to 3,072 samples per week**



Extracted Viral RNA Samples

cDNA Synthesis & Probe Hybridization

Circularization & Cleanup

PCR with Asymmetric Barcoded M13 Primers (F/R)

Pool Barcoded Samples

SMRTbell Library Construction & Sequencing Preparation

SMRT Sequencing (Sequel II or IIe System)



HiFiViral SARS-CoV-2 Data Analysis With SMRT Link

4 – 16 h

4 – 5 h

6 h

14 – 16 h**




* HiFiViral SARS-CoV-2 Kit demonstrated use cases include RNA-extracted samples such as nasopharyngeal or saliva swabs from human SARS-CoV-2 PCR+ cohort samples.

** For multi-SMRT Cell runs, sequencing + data analysis time is ~14 – 16 h for the first cell. For subsequent SMRT Cells, sequencing + data analysis time is ~9 – 10 h per cell

HiFiViral SARS-CoV-2 KIT USES *MOLECULAR INVERSION PROBE* TECHNOLOGY FOR EFFICIENT VIRAL GENOME ENRICHMENT


Overview of MIP-Based Viral Enrichment Enzymatic Reaction Steps

- Single pool of 969 ssDNA MIPs comprised of two probe arms connected by a common linker (30 bp)
- MIPs tile SARS-CoV-2 genome at 22-fold target coverage
- Capture of 675-bp target sequences is performed by circularization of MIPs via "fill-in" enzymatic reaction
- PCR amplification using universal (M13) primers adds unique (dual index) asymmetric barcodes to each sample to enable multiplexed analyses




HiFiViral SARS-CoV-2 KIT USES *MOLECULAR INVERSION PROBE* TECHNOLOGY FOR EFFICIENT VIRAL GENOME ENRICHMENT (CONT.)

Dense MIP-Based Tiling of Target Sequences Enables Robust Coverage



↓ PCR with M13 universal primers to add dual indices for sample multiplexing




Advantages of MIPs

- **Higher Specificity**
 - Each MIP molecule contains two probe arms
- **Easier Workflow**
 - Unlike traditional PCR-based targeting with overlapping primers, overlapping MIPs can be used in a single reaction leading to fewer plates and fewer touch points
- **More Robust Probe Design**
 - ~1000 ssDNA probes tile target SARS-CoV-2 genome at 22-fold coverage
 - More tolerant to viral RNA sample degradation and a wider range of input RNA quantities
 - Resilient to mutation-induced probe dropouts with new viral genomic variants

Asymmetrically Barcoded Double-Stranded Library Molecules (~800 bp)*

* Not to scale.

HIFIVIRAL SARS-COV-2 SEQUENCING REQUIRES FEWER READS FOR COMPLETE VIRAL GENOME COVERAGE




TECHNOLOGY	# OF READS FOR COMPLETE COVERAGE	MINIMUM READ DEPTH
PacBio HiFi	1,000	4-fold
Oxford Nanopore	10,000	20-fold
Illumina	1,000,000	10-fold

**HiFi reads are more accurate
→ Fewer reads simplifies analysis**

96-Plex of Twist Control Samples

Samples with Ct ≤ 30 achieved complete genome coverage* with 1000 HiFi reads



High Viral Copy Number Abundance

Low Viral Copy Number Abundance

* Complete = ≥90% genome coverage



HiFiViral SARS-CoV-2 Kit Workflow Overview

HIFIVIRAL SARS-COV-2 KIT LIBRARY PREPARATION PROCEDURE DESCRIPTION

- **Procedure & Checklist – HiFiViral for SARS-CoV-2 Workflow** ([102-188-800](#)) describes a viral enrichment and library preparation procedure for whole viral genome sequencing of multiplexed SARS-CoV-2 samples on the Sequel II and IIe Systems using HiFiViral SARS-CoV-2 Kit ([102-132-000](#)) and SMRTbell Template Prep Kit 2.0 ([100-938-900](#))




HiFiViral SARS-CoV-2 Kit
([102-132-000](#))



SMRTbell Express TPK 2.0
([100-938-900](#))

- This procedure utilizes **molecular inversion probe** (MIP)-based chemistry to enrich the SARS-CoV-2 genome with tiled probes that create highly-redundant overlapping amplicons, which are barcoded and pooled for construction into a single SMRTbell library for sequencing
- Viral enrichment uses an **addition-only 4-step workflow with color-coded master mixes** to simplify setup
- End-to-end workflow from cDNA synthesis through to SMRTbell library construction, sequencing & analysis can be completed in as short as 28 – 42 hours depending on desired hybridization time



PacBio Documentation ([102-188-800](#))

RESEARCH FOCUS
MICROBIOLOGY AND INFECTIOUS DISEASE
PacBio COVID-19 Sequencing Tools and Resources



HiFiViral SARS-CoV-2 KIT PRODUCT DESCRIPTION

HiFiViral SARS-CoV-2 Kit (102-132-000)

- Assay kit designed for targeted enrichment and barcoding of up to 384 human SARS-CoV-2-positive samples for full-length viral genomic sequencing on PacBio Sequel II or Ilx Systems
- Kit contains two components: 1) SARS-CoV-2 Enrichment Kit; and 2) Barcoded M13 Primer Plate

SARS-CoV-2 Enrichment Kit

- The SARS-CoV-2 Enrichment Kit contains all reagents for enrichment using Molecular Inversion Probes (MIPs) of extracted RNA virus from cohort samples infected with the SARS-CoV-2 virus. This kit is to be used in conjunction with the Barcoded M13 Primer Plate.
- The results of the kit are enriched DNA fragments of ~800 bp in length that can be used to prepare a SMRTbell library for sequencing.
- Reagent quantities support preparation of 384 samples with flexible scaling down to batches of 24 samples.

Barcoded M13 Primer Plate*


- 1 premixed primer plate containing 384 barcoded M13 primer pairs for asymmetric (dual index) barcoding of multiplexed SMRTbell libraries
- Single-use per well with pierceable foil (can reseal between sample batches)



HiFiViral SARS-CoV-2 Kit (102-132-000)






SARS-CoV-2 Enrichment Kit



Barcoded M13 Primer Plate

* Barcoded M13 Primer Plate (102-135-500) part may also be ordered separately for use with other multiplexed SMRT Sequencing applications.

HiFiViral SARS-CoV-2 KIT COMPONENTS

KIT PRODUCT OR COMPONENT	KIT SUBCOMPONENT	PART NUMBER	QUANTITY	NO. OF REACTIONS SUPPORTED
<p>HiFiViral SARS-CoV-2 Kit (NEW)</p> 		102-132-000		384
<p>SARS-CoV-2 Enrichment Kit</p> 	<ul style="list-style-type: none"> 1 Probe Mix 2 Fill-In Mix 3 Cleanup Mix 4 Reverse Transcriptase Mix 5 PCR Mix 		<ul style="list-style-type: none"> 1 Tube 1 Tube 1 Tube 1 Tube 3 Tubes 	
<p>Barcoded M13 Primer Plate*</p> 	Premixed Primer Plate		1 Plate	

* Barcoded M13 Primer Plate ([102-135-500](#)) part may also be ordered separately for use with other multiplexed SMRT Sequencing applications.

HiFiViral SARS-CoV-2 KIT WORKFLOW OVERVIEW

Sample RNA Extraction


- Isolate viral RNA from human SARS-CoV-2 PCR+ samples using third-party protocols



Sample RNA Extraction

Viral Genome Enrichment with HiFiViral SARS-CoV-2 Kit

- Addition-only workflow features a visible color change with each reagent addition step to signal success
- All reactions performed on one sample plate
- MIP-based viral enrichment workflow times:
 - Overnight cDNA Synthesis (1 h) + Hybridization (16 h) with option to reduce hybridization time to 4 h for faster turnaround time. (A longer hybridization time boosts HiFi Data Yield for high-Ct samples)
 - MIP Circularization, Cleanup, PCR and Pooling steps can be completed in ~4 – 5 h
- Amplify and asymmetrically barcode up to 384 SARS-CoV-2 samples (per SMRT Cell 8M) for multiplexing in a single library using PacBio-Barcoded M13 Primers



cDNA Synthesis & Probe Hybridization

Circularization & Cleanup

PCR with Asymmetric Barcoded M13 Primers (F/R)


Pool Barcoded Samples

4 – 16 h

4 – 5 h

SMRTbell Library Construction & Sequencing Prep

- Library Prep: SMRTbell Express TPK 2.0; SMRTbell Enzyme Cleanup Kit 2.0
- Sequencing Prep: Sequencing Primer v5; Binding Kit 2.1; ProNex Bead Cleanup



SMRTbell Library Construction & Sequencing Preparation

6 h

Sequencing & Data Analysis

- Use 8-h movie collection time per SMRT Cell 8M
- Load up to 8 SMRT Cells per run to analyze up to 3,072 samples in 1 week
- Use SMRT Link HiFiViral SARS-CoV-2 analysis application for data analysis



SMRT Sequencing & Analysis

14 – 16 h



Multiplexed Library Preparation Using Molecular Inversion Probe-Based Enrichment with the HiFiViral SARS-CoV-2 Kit

PROCEDURE & CHECKLIST – PACBIO HIFIVIRAL HIGH-THROUGHPUT MULTIPLEXING FOR FULL-VIRAL GENOME SEQUENCING OF SARS-COV-2

Procedure & Checklist [102-188-800](#) describes a viral enrichment and library preparation procedure for whole viral genome sequencing of multiplexed SARS-CoV-2 samples on the Sequel II and IIe Systems using HiFiViral SARS-CoV-2 Kit (102-132-000) and SMRTbell Template Prep Kit 2.0 (100-938-900)




HiFiViral SARS-CoV-2 Kit
([102-132-000](#))



SMRTbell Express TPK 2.0
([100-938-900](#))

Protocol Contents

1. RNA input requirements and best practices recommendations for preparing master mixes, handling RNA samples, and sealing reaction plates
2. Instructions for performing enrichment of SARS-CoV-2 viral cDNA products using HiFiViral SARS-CoV-2 Kit ([102-132-000](#))
3. Instructions for pooling amplified SARS-CoV-2 cDNA products and constructing SMRTbell libraries using SMRTbell Express Template Prep Kit 2.0 ([100-938-900](#))



REQUIRED MATERIALS & EQUIPMENT

ITEM	WHERE USED	VENDOR	PART NUMBER
RNA Preparation			
Nuclease-Free Water	RNA Preparation	Any	Vendor-specific
RNaseZap	RNA Preparation	Thermo Fisher Scientific	AM9780
SARS-CoV-2 RNA Viral Enrichment			
HiFiViral SARS-CoV-2 Kit (Includes items below) SARS-CoV-2 Enrichment Kit Barcoded M13 Primer Plate	cDNA Synthesis & Probe Hybridization Reaction Fill Reaction Cleanup Reaction PCR & Barcoding Reaction Library Pooling	PacBio	102-132-000
SMRTbell Library Preparation			
SMRTbell Express Template Prep Kit 2.0	Library Construction	PacBio	100-938-900
SMRTbell Enzyme Cleanup Kit 2.0	Library Construction	PacBio	101-932-600
DynaMag-2 Magnet	Library Purification	Thermo Fisher Scientific	12321D
Absolute Ethanol, Molecular Biology or ACS Grade	Library Purification	Any	Vendor-specific
ProNex Beads	Library Purification	Promega	NG2001-10mL / NG2002-125mL / NG2003-500mL
DNA LoBind Tubes	Library Construction	Eppendorf	22431021 (1.5 mL) / 22431048 (2.0 mL)

REQUIRED MATERIALS & EQUIPMENT (CONT.)

ITEM	WHERE USED	VENDOR	PART NUMBER
Sequencing Preparation			
Sequel II Binding Kit 2.1 and Int Ctrl 1.0	Sequencing on the Sequel II and IIe Systems	PacBio	101-843-000
Sequel II Sequencing Kit 2.0	Sequencing on the Sequel II and IIe Systems	PacBio	101-820-200
SMRT Cell 8M Tray	Sequencing on the Sequel II and IIe Systems	PacBio	101-389-001
Sequel SMRT Oil	Sequencing on the Sequel II and IIe Systems	PacBio	100-621-300
Sequel Pipette Tips v2	Sequencing on the Sequel II and IIe Systems	PacBio	100-667-601
Sequel Mixing Plates	Sequencing on the Sequel II and IIe Systems	PacBio	100-667-500
Sample Plate	Sequencing on the Sequel II and IIe Systems	PacBio	HSP9601
Tube Septa	Sequencing on the Sequel II and IIe Systems	PacBio	001-292-541
Sequel Sample Plate Foil	Sequencing on the Sequel II and IIe Systems	PacBio	100-667-400
DNA / Library QC Evaluation			
Qubit 4 Fluorometer	DNA Quantitation	Thermo Fisher Scientific	Q33238
Qubit 1x dsDNA HS Assay Kit	DNA Quantitation	Thermo Fisher Scientific	Q33230
Bioanalyzer 2100	DNA Sizing	Agilent	G2939A
Agilent DNA 12000 Kit	DNA Sizing	Agilent	5067-1508

REQUIRED MATERIALS & EQUIPMENT (CONT.)

ITEM	WHERE USED	VENDOR	PART NUMBER
General			
96-well PCR plates	SARS-CoV-2 RNA Viral Enrichment	Bio-Rad	HSP9601
Microseal 'B' Film	SARS-CoV-2 RNA Viral Enrichment	Bio-Rad	MSB1001
Film sealing roller for PCR plates	SARS-CoV-2 RNA Viral Enrichment	Bio-Rad	MSR0001
Thermal Cycler With Heated Lid (Examples below) VeritiPro Thermal Cycler, 96 well ProFlex PCR System	SARS-CoV-2 RNA Viral Enrichment	Thermo Fisher Scientific	A48141 (VentiPro) / 4483636 (ProFlex)
PCR Tube Strips, 0.2 mL	SARS-CoV-2 RNA Viral Enrichment	USA Scientific	1402-4708
96-Well Plate Centrifuge	SARS-CoV-2 RNA Viral Enrichment	Any Vendor	Vendor-specific
8- or 12-Multichannel Pipette	SARS-CoV-2 RNA Viral Enrichment	Any Vendor	Vendor-specific

HiFiViral SARS-CoV-2 KIT WORKFLOW DETAILS



1 Viral Genome Enrichment with HiFiViral SARS-CoV-2 Kit (102-132-000)

i. cDNA Synthesis & Hybridization With Molecular Inversion Probes (MIPs)

- For input into cDNA synthesis & probe hybridization, aim to start with $\geq 10,000$ copies of RNA
- For SARS-CoV-2 viral enrichment using ssDNA MIPs, a 16-h hybridization time (55°C) is recommended for high-Ct samples (Ct >25). Optionally, can reduce hybridization time to as short as 4 h for low-Ct samples



ii. Circularization of MIPs

- Circularize MIPs bound to target SARS-CoV-2 sequences via fill-in reaction between ssDNA probe arms



iii. Enzymatic Cleanup

- Remove unbound MIPs via enzymatic cleanup




iv. PCR Amplification with Asymmetric Barcoded M13 Primers

- Amplify captured SARS-CoV-2 target sequences via 1-step PCR with asymmetric barcoded M13 primers (included with HiFiViral SARS-CoV-2 Kit) to produce barcoded ~800 bp dsDNA amplicon products



v. Pool Barcoded Samples

- Perform equal-volume pooling of barcoded SARS-CoV-2 amplicon products from up to 384 cohort samples (per SMRT Cell 8M) and purify the resulting mixture using ProNex Beads.
- After ProNex Bead cleanup, the pooled sample is used for construction into a single, multiplexed SMRTbell library



HiFiViral SARS-CoV-2 KIT WORKFLOW DETAILS (CONT.)



2 SMRTbell Library Construction with SMRTbell Express Template Prep Kit 2.0 (100-938-900)

- The amount of total pooled (barcoded) DNA required for SMRTbell library construction is 500 ng – 1000 ng.
- Library construction can be completed in ~5 h
- Typical library construction yield is ≥40%




3 SMRT Sequencing

- Anneal Sequencing Primer v5 (15 min), bind Sequel Polymerase 2.2 (15 min), and perform complex cleanup with ProNex Beads (0.5 h)
- Use 8-h movie collection time per SMRT Cell 8M for sequencing HiFiViral SARS-CoV-2 Kit samples on the Sequel II or IIe System
- Load up to 8 SMRT Cells per run to analyze up to 3,072 samples in 1 week



4 Data Analysis

- Use SMRT Link HiFiViral SARS-CoV-2 analysis application to perform variant calling




* For multi-SMRT Cell runs, sequencing + data analysis time is ~14 – 16 h for the first cell. For subsequent SMRT Cells, sequencing + data analysis time is reduced to ~9 – 10 h per cell due to parallelization of sequencing and analysis functions during the instrument run.

HiFiViral SARS-CoV-2 KIT SAMPLE PREP WORKFLOW TIMING SUMMARY

Efficient Workflow (~1 – 2.5 Hours Hands-On Time) Enables Sample to Answer in ~28 – 42 Hours

WORKFLOW STEP	HANDS-ON (MIN)	WALK-AWAY (HRS)
SARS-CoV-2 RNA Enrichment (~22 h)		
cDNA Synthesis	5 – 15	1.0
Probe Hybridization with MIPs	5 – 15	4.0 – 16.0
Circularization (Fill-in Reaction)	5 – 15	1.0
Enzymatic Cleanup Reaction	5 – 15	1.2
PCR with Barcoded M13 Primers	10 – 30	1.5
Pooling (DNA sizing QC is optional)	5 – 10	—
1.3X ProNex Bead Cleanup + Qubit Assay	5 – 10	0.3
Total	~40 – 110	~9.0 – 21.0
SMRTbell Library Construction (~5 h)		
DNA Damage Repair	2 – 4	0.5
End Repair / A-Tailing	2 – 4	1.0
Adapter Ligation	2 – 4	1.2
1.3X ProNex Bead Cleanup	2 – 4	0.3
Nuclease Treatment	2 – 4	0.5
1.3X ProNex Bead Cleanup + Qubit Assay	5 – 10	0.3
Total	~15 – 30	~3.8
Sequencing Preparation (~1.5 h)		
Anneal Sequencing Primer	2.5 – 5	0.25
Bind Polymerase	2.5 – 5	0.25
1.2X ProNex Bead Complex Cleanup	5 – 10	0.5
Total	~10 – 20	~1.0



* For multi-SMRT Cell runs, sequencing + data analysis time is ~14 – 16 h for the first cell. For subsequent SMRT Cells, sequencing + data analysis time is reduced to ~9 – 10 h per cell due to parallelization of sequencing and analysis functions during the instrument run.

RNA INPUT REQUIREMENTS FOR VIRAL ENRICHMENT USING HiFiViral SARS-CoV-2 KIT

- Best results will be achieved if reactions contain **at least 10,000 copies of RNA**.
 - Samples with higher copy numbers of RNA virus will generally produce superior results.
 - See at table at right for example viral copy number values converted from a Ct scale*
- Purified RNA should be resuspended in RNase-free water or TE with a pH no greater than 7.5.
- Contaminants including ethanol, sodium azide, sodium acetate, and guanidine salts may affect performance.
- DNase treatment is optional but the presence of small amounts of human DNA should not affect performance.
- If RNA is quantified, a method that is specific for RNA is recommended (e.g., Qubit RNA BR Assay Kit or qRT-PCR), rather than one that will also detect DNA.
- To reduce inter-sample performance variability, all samples in a batch should be quantified using the same method and normalized to the same concentration.

Example viral copy number values shown in Table below are converted from a Ct scale after Han *et al.* 2021.

Sample Ct	Viral Copy Number*
19	6 Million
20	3 Million
21	1 Million
24	100,000
27	10,000
30	1,000
33	100
35	3

* **NOTE:** A Ct value itself **cannot** be directly interpreted as viral load without a standard curve using reference materials. [See Han M.S., et al. (2021). RT-PCR for SARS-CoV-2: quantitative versus qualitative. *The Lancet Infectious Disease* 21(2) p165]

GENERAL BEST PRACTICES RECOMMENDATIONS FOR VIRAL ENRICHMENT USING HiFiViral SARS-CoV-2 KIT

Best Practices

Master Mixes

1. Prepare master mixes in a PCR workstation.
2. The PCR workstation should be UV-irradiated after each setup. If unsure, UV-irradiate the workstation before setting up a master mix.
NOTE: do not turn on the UV light when reagents are in the workstation.
3. Master mixes are prepared in 0.5mL, 1.5 mL or 2 mL PCR tubes.
down.
4. If using multichannel pipette to transfer master mixes, pre-aliquot appropriate volume with overage into PCR strip tubes.

Samples

1. RNA samples should be stored at -80°C until use and thawed on ice.
2. Heavily degraded RNA or RNA samples with many freeze-thaw cycles should be avoided.
3. All work surfaces and gloves should be sanitized with RNaseZap (or the equivalent) prior to setting up.
4. For most consistent performance, all samples included in a batch should be from the same sample type and extracted by the same RNA extraction procedure.
5. A no-RNA control is recommended but not required.
6. Upon thawing frozen samples, briefly vortex and spin down prior to use.

Reaction Plates

1. Always seal plates with Microseal 'B' Film (clear adhesive). Foil seals are not recommended for any step in this protocol. However, they can be used for plates that will be placed in the freezer for storage.
2. Using a roller for Microseal 'B' Film, apply firm pressure and seal over the tops of all wells. Ensure all wells, especially those along the edges of the plate, are visibly sealed.
3. Inspect the corners of the plate to confirm that the seal is in contact with the plate. If not, apply firm pressure and roll until the film is in contact with the plate.
4. When removing plate seals, a heated plate sealer is recommended to loosen the adhesive.
5. Centrifuge in an Eppendorf 5810 fitted with a swinging bucket plate rotor at maximum rpm for approximately 30 sec.
6. After centrifugation, inspect the bottom of the plate to ensure the expected volume is present in every well.

Preparing Master Mixes

- Prepare master mixes in a PCR workstation and have them ready before the end of the prior incubation steps
- Use multichannel or electronic pipettes to facilitate transfer of master mixes to sample wells



Handling RNA Samples

- Use special care when handling small volumes of reagents
- Be careful when removing plate seals to avoid cross contamination

Preparing Reaction Plates

- Always perform a visual check of liquid volumes before and after each incubation step
- Verify that the liquid solution color at each reach step is correct
- Proper plate sealing is **critical**, especially for the overnight probe hybridization step

HiFiViral SARS-CoV-2 KIT PROCEDURAL NOTES

1. cDNA Synthesis and Probe Hybridization

cDNA Synthesis and Probe Hybridization

Before setting up the reaction, the workstation should be sanitized with RNaseZap and UV-irradiated without the presence of the reagents. All samples and reagents should be kept on ice while setting up the reaction.

- Prepare labware and reagents.
 - Label one or more 96-well PCR plates. Alternatively, for a small number of reactions, PCR tube strips may be used.
 - Retrieve extracted RNA samples from storage.
- Pipette 6 μL of sample RNA into each well of the reaction plate. Be sure to follow RNA input volume recommendations. Use nuclease-free water to adjust sample RNA volume, if needed.
- Prepare RT-Hybridization Master Mix on ice.
 - Allow RT Mix and Probe Mix to fully thaw. Briefly vortex and spin down.
 - Prepare master mix with 12.5% overage as indicated in the table below. Preparing fewer than 24 reactions at a time is not recommended.
 - RT Mix is viscous, pipette slowly.

Reagent	1X reaction	96 reactions (with overage)	✓	Notes
RT Mix	1.6 μL	172.8 μL		
Probe Mix	0.4 μL	43.2 μL		

- Transfer 2 μL of RT-Hybridization Master Mix into each sample-containing well in the reaction plate. RT-Hybridization Master Mix is viscous, pipette slowly.
- Seal the plate tightly with a film. Poor sealing could result in significant evaporation.
- Spin down the 96-well plate.
- Vortex a few times with short pulses and spin down.
- Perform a quick visual check of the liquid level and take note of any well with low volume. The reaction should now be a homogenous **pale blue** color.
- Place the reaction plate in the thermal cycler and run the following program (set the heated lid at 105°C).

Step	Temperature	Time
1	25°C	10 minutes
2	50°C	50 minutes
3	95°C	
4	55°C	24 hours * (16hrs for Probe Hybridization and 1hr for Fill-in Reaction)
5	55°C	Hold

- Make a note of the thermal cycler start time. A hybridization time of 16 hours (the 55°C step) is recommended for high Ct samples (Ct >25). A 4hr hybridization could be considered if most of samples have low Ct value (Ct <25). Start preparing for the fill reaction just prior to the end of hybridization (approximately 17 hours from the start of the cycling program).


Preparing Master Mixes

- Slowly pipette small reaction volumes and viscous reagents (Master Mix volumes shown in the table only include 12.5% overage)



Preparing Reaction Plates

- Seal reaction plates tightly with Microseal 'B' Film to minimize evaporation, especially along the plate edges and corners
- Verify that the liquid solution color for each Hybridization reaction is **blue** and homogeneous



Starting and Monitoring Hyb Reactions

- Do not use questionable or problematic thermal cycler equipment for this viral enrichment workflow
- A 16-hour hybridization time is recommended – Make note of the reaction start time (incubating slightly longer than 16 hours should not have a negative impact)
- Keep the thermal cycler program running after probe hybridization is completed to maintain proper temperature control of the heating block

HiFiViral SARS-CoV-2 KIT PROCEDURAL NOTES (CONT.)

2. Circularization (Fill Reaction)

Fill Reaction

Before the end of the probe hybridization reaction, allow the Fill-in Mix to fully thaw. Briefly vortex and spin down. Do not remove the reaction plate from the thermal cycler until the reagent is ready and the hybridization time is over. Correct timing is important to maximize result quality.

1. Remove the sample plate from the thermal cycler. **Keep the program running.**
2. Spin down the plate, perform a quick visual check of the liquid level to make sure there are no droplets on the top seal or side walls, and remove the seal carefully to avoid cross contamination.
3. At room temperature, transfer **2 μ L** of Fill-in Mix to each sample well.
It is important to perform the transfer as fast as possible to minimize non-specific binding; aim to finish within 5 minutes.
4. Reseal the plate tightly with a new film, vortex a few times with short pulses, and spin down the plate.
5. Perform a quick visual check of the liquid level and take note of any well with low volume. The reaction should now be a homogenous **pale green** color.
6. Place the reaction plate in the thermal cycler and continue the program for another 60 minutes.
The time the reaction plate was returned to the thermal cycler; correct timing is important to maximize result quality.

Preparing Fill Reaction Plates

- Add reagents at room temperature, **DO NOT** cool on ice
- Fill Reaction steps are time sensitive – Work quickly with a multichannel pipettor to complete all liquid transfer steps within 5 minutes for best capture results
- Verify that the liquid solution color for each Fill Reaction is **green** and homogeneous



3. Cleanup Reaction

Cleanup Reaction

Before the end of the fill reaction, allow the Cleanup Mix to fully thaw. Briefly vortex and spin down. Do not remove the reaction plate from the thermal cycler until the reagent is ready. Correct timing is important to maximize result quality.



1. Remove the sample plate from the thermal cycler.
2. Spin down the plate, perform a quick visual check of the liquid level to make sure there are no droplets on the top seal or side walls, and remove the seal carefully to avoid cross contamination.
3. At room temperature, transfer **2 μ L** of Cleanup Mix to each sample well.
It is important to perform the transfer to minimize non-specific binding; aim to finish within 10 minutes.
4. Reseal the plate tightly with a new film, vortex a few times with short pulses, and spin down the plate.
5. Perform a quick visual check of the liquid level and take note of any well with low volume. The reaction should now be a homogenous **red** color.
6. Place the reaction plate in the thermal cycler and run the following program (set the heated lid at 105°C).

Step	Temperature	Time
1	45°C	60 minutes
2	95°C	3 minutes
3	4°C	Hold

7. The program is approximately 65 minutes to run; proceed immediately to the cDNA amplification step.

Preparing Cleanup Reaction Plates

- Add reagents at room temperature, **DO NOT** cool on ice
- Cleanup Reaction steps are time sensitive – Work quickly with a multichannel pipettor to complete all liquid transfer steps within 5 minutes for best capture results
- Verify that the liquid solution color for each Cleanup Reaction is **red** and homogeneous



HiFiViral SARS-CoV-2 KIT PROCEDURAL NOTES (CONT.)

4. PCR Amplification With Barcoded M13 Primers

cDNA Amplification

Before the end of the cleanup reaction, allow the PCR Mix and the primer plate to fully thaw. Spin down the primer plate before opening. Briefly vortex the PCR Mix and spin down. The reaction plate and reagents should be kept on ice while setting up the reaction.

1. Remove the sample plate from the thermal cycler.
2. Spin down the plate, perform a quick visual check of the liquid level to make sure there are no droplets on the top seal or side walls, and remove the seal carefully to avoid cross contamination.
3. Prepare the PCR reaction as follow.

Reagent	Stock Conc.	1X reaction	
Cleanup reaction		9.6 μ L	
PCR Mix		12 μ L	
Asymmetric Barcoded M13 Primer Mix	10 μ M	2.4 μ L	
Total Volume		24 μ L	

* The expected volume after the cleanup reaction is approximately 9.6 μ L, considering some degree of evaporation during the prior steps


4. Using a multichannel pipette, transfer 12 μ L of PCR Mix to the sample plate.
5. Transfer 2.4 μ L of premixed asymmetric barcoded M13 primers from the primer plate to the corresponding sample wells.
6. The total reaction volume in each well is approximately 24.0 μ L.
7. Reseal the plate tightly with a new film, vortex a few times with short pulses, and spin down the plate.
8. Perform a quick visual check of liquid level and take note of any well with low volume. The reaction should now be a homogenous **magenta** color.
9. Place the PCR reactions in a thermal cycler and run the following program (set the heated lid at 105°C).

Step	Temperature	Time
1	95°C	3 minutes
2	98°C	15 seconds
3	55°C	30 seconds
4	72°C	1 minute 30 seconds
5	Repeat steps 2 to 4 (26 times)	
6	4°C	Hold

10. After amplification, briefly spin down the plate.
11. Immediately proceed to the "Sample Pooling for Library Construction" section if not performing the optional Library Quantitation/QC step. Alternatively, the reaction plate can be stored at -20°C until further processing.

Preparing PCR Reaction Plates

- Expected sample volume after cleanup step is ~9.6 μ L.
- PCR amplification step is not time-sensitive
- Verify that the liquid solution color for each PCR Reaction is **magenta** and homogeneous



Starting and Monitoring PCR Reactions


- PCR thermal cycler program at this step takes ~1.5 hours to complete (27 cycles)
- Expect some degree of cumulative evaporation loss to occur from completing previous steps in the workflow – If any sample in a well has significantly less than 9.6 μ L, add nuclease-free water to bring up the sample volume and document this action
- After completing the PCR step, amplified cDNA samples can be stored at -20°C until further processing

HiFiViral SARS-CoV-2 KIT PROCEDURAL NOTES (CONT.)

4. PCR Amplification With Barcoded M13 Primers (Cont.)

Library Quantification/QC (Optional)

1. Remove the reaction plate from the thermal cycler.
2. Spin down the reaction plate and perform a quick visual check of the liquid level. Take note of any well with low volume, which indicates excessive evaporation during amplification.
3. Remove the seal carefully to avoid cross contamination.
4. Use 1 μ L of sample to quantify with a Qubit dsDNA HS kit.
5. Individual sample QC can be performed on the Agilent 2100 Bioanalyzer. Use a DNA12000 chip and follow the manufacturer's setup instruction.
6. A target peak of ≥ 700 bp should be detected. A small peak of ~ 170 -200 bp representing non-specific amplicons may or may not be present. The ~ 170 -200 bp amplicons will be removed when the sample pool is purified.




Example post-PCR DNA sizing analysis results for extracted viral RNA samples.

- Spot-checking PCR amplification products prior to pooling is highly recommended when performing the HiFiViral workflow for the first time
- 1.3X ProNex Bead purification can help remove non-specific amplification products

Post-PCR DNA Quantification and DNA Sizing QC (Optional)


Performing **post-PCR DNA sizing quantification and sizing QC steps is recommended** and can be useful for verifying sample integrity prior to SMRTbell library construction as well as downstream troubleshooting



Example post-PCR DNA sizing analysis results for ATCC Control RNA samples.

- Going from high to low copy number, the on-target band diminishes, and the amount of non-specific amplification products increases
- 1.3X ProNex Bead purification can help remove non-specific amplification products

PCR Product Yield vs. Input Control RNA Copy Number



Example post-PCR yield results for ATCC Control RNA samples.

- Higher-copy number samples are generally correlated with higher PCR yields (*via* Qubit dsDNA HS assay quantification)

HiFiViral SARS-CoV-2 KIT PROCEDURAL NOTES (CONT.)


5. Sample Pooling for SMRTbell Library Construction

Sample Pooling for Library Construction

1. Remove the reaction plate from the thermal cycler.
2. Spin down the reaction plate and perform a quick visual check of the liquid level. Take note of any well with low volume, which indicates excessive evaporation during amplification.
3. Remove the seal carefully to avoid cross contamination.
4. Transfer a minimum of 5 μL per reaction into a clean 1.5 mL or 2.0 mL Lo-bind tube.
 - a. The total pool should be at least 100 μL . For example, if running 8 reactions, pool 12.5 μL per reaction.
 - b. If pooling 384 reactions, vortex to mix and transfer no more than 800 μL to a new 2.0 mL Lo-bind tube for purification. Save the rest of the sample pool at -20°C .

Preparing Samples for Pooling

- Transfer a minimum of 5 μL per reaction into a clean 1.5 mL or 2.0 mL Lo-bind tube.
- The total pool volume should be at least 100 μL
 - If running 8 samples, pool 12.5 μL from each PCR reaction
 - For 96 samples, Total Pool Volume = 480 μL
 - For 384 samples, Total Pool Volume = 1920 μL
- Note: If pooling 384 reactions, the total volume is too large for a 1.5 mL tube
 - Transfer no more than 800 μL to a new 2.0 mL Lo-bind tube for purification. (Save the rest of the sample pool at -20°C .)



HiFiViral SARS-CoV-2 KIT PROCEDURAL NOTES (CONT.)

5. Sample Pooling for SMRTbell Library Construction (Cont.)

Purification of Pooled Library

STEP	✓	Purification with ProNex Beads	Notes
1		Add 1.3X volume of resuspended, room-temperature ProNex beads to the pooled library. Pipette mix 10 times. Perform a quick-spin to collect all liquid from the sides of the tube.	
2		Incubate the sample on the bench top for 5 minutes at room temperature.	
3		Place the tube on a magnetic stand to separate the beads from the supernatant. Use a P200 pipette to remove the supernatant.	
4		Wash 2 times with 1400 μ L (or enough to fully cover the beads) of freshly prepared 80% ethanol. After removal of the second 1400 μ L ethanol wash, spin the tube and return to the magnetic stand, and remove residual ethanol with a P200 pipette.	
5		Remove the tube from the magnetic stand. Immediately pipette 100 μ L of the supernatant and pipette the mix to resuspend. Perform a quick-spin to collect all liquid from the sides of the tube. Place at room temperature for 5 minutes to elute the DNA from the beads.	
6		Place the tube on a magnetic stand to separate the beads from the supernatant. Transfer the eluted DNA sample to a new tube.	
7		Use 1 μ L of sample to quantify with a Qubit dsDNA HS kit.	

Purifying Pooled Samples



- Add 1.3X volume of resuspended, room-temperature ProNex beads to the pooled library.
 - Bead incubation: 5 mins, Room Temperature
 - Elution incubation: 5 mins, Room Temperature
- The total amount of purified pooled (barcoded) DNA required for SMRTbell library construction is 500-1000 ng.

Library Quantification/QC (Optional)

1. Pooled sample QC can be performed on the Agilent 2100 Bioanalyzer. Use a DNA12000 chip and follow the manufacturer's setup instruction.
2. A target peak of ≥ 700 bp should be detected. Non-specific amplicons (~170-200 bp) should be removed completely.

DNA Sizing QC

- DNA sizing QC can optionally be performed on the pooled sample using an Agilent 2100 Bioanalyzer
 - A target peak of ≥ 700 bp should be detected
 - Non-specific amplicons (~170-200 bp) should be removed completely.



SMRTBELL EXPRESS TEMPLATE PREP KIT 2.0 AND SMRTBELL ENZYME CLEANUP KIT 2.0 REAGENT HANDLING RECOMMENDATIONS

- Several reagents in the kit are sensitive to temperature and vortexing
- PacBio highly recommends:
 - Never leaving reagents at room temperature
 - Working on ice at all times when preparing master mixes for SMRTbell library construction
 - Finger tapping followed by a quick-spin prior to use

SMRTbell Express TPK 2.0
([100-938-900](https://www.pacb.com/products-services/sequencing-chemistry/sequencing-chemistry-products/sequencing-chemistry-products-smrtbell-express-tpk-2-0/))



LIST OF TEMPERATURE-SENSITIVE REAGENTS INCLUDED IN SMRTBELL EXPRESS TPK 2.0 AND SMRTBELL ENZYME CLEANUP KIT 2.0.

PACBIO KIT	REAGENT	WHERE USED
SMRTbell Express Template Prep Kit 2.0 (PN 100-938-900)	DNA Prep Additive	Remove Single-Strand Overhangs
	DNA Prep Enzyme	Remove Single-Strand Overhangs
	DNA Damage Repair Mix v2	DNA Damage Repair
	End Prep Mix	End-Repair/A-tailing
	Overhang Adapter v3	Ligation
	Ligation Mix	Ligation
	Ligation Additive	Ligation
Ligation Enhancer	Ligation	
SMRTbell Enzyme Cleanup Kit 2.0 (PN 101-932-600)	SMRTbell Enzyme Clean Up Mix	Nuclease Treatment
	SMRTbell Enzyme Clean Up Buffer	Nuclease Treatment



Multiplexed SARS-CoV-2 Library Sequencing Workflow Recommendations



SAMPLE SETUP AND RUN DESIGN RECOMMENDATIONS FOR HiFiViral SARS-CoV-2 LIBRARY SAMPLES (SEQUEL II/IE SYSTEMS)

Follow **SMRT Link Sample Setup** instructions using the recommendations provided in **Quick Reference Card – Loading and Pre-Extension Time Recommendations for the Sequel II/IE Systems** ([101-769-100](tel:101-769-100)) for preparing HiFiViral samples for sequencing

Application	Library Prep Kit*	Sequencing Primer (Annealing Time)	Binding Kit (Binding Time)	Complex Cleanup	Loading Concentration Range (pM)
HiFiViral SARS-CoV-2 (1 kb)	Express TPK 2.0	v4 (15 min)	Binding Kit 2.1 (15 min)	1.2X ProNex Beads	100 - 300

Application	Pre-Extension Time (hours)	Adaptive Loading Target (P1 + P2)	Movie Collection Time (hours)
HiFiViral SARS-CoV-2 (1 kb)	0	N/A	8

→ For **SMRT Link v10.2**: Select **'Viral Sequencing / HiFiViral SARS-CoV02'** from the **Application** field drop-down menu in the SMRT Link Sample Setup and SMRT Link Run Design user interface



IMPORTING THE BARCODE FASTA FILE INTO SMRT LINK FOR AUTOMATED DEMULTIPLEXING OF HiFiViral SARS-CoV-2 LIBRARY SAMPLES

Note: SMRT Link v10.2 software installations by default come **pre-bundled** with FASTA files containing a list of PacBio barcodes recommended for use with specific multiplexed SMRT sequencing applications

If your SMRT Link installation does **not** already include an appropriate barcode FASTA file, the following steps describe how to import such a file for use in automated demultiplexing (refer to “Importing Data” section in the [SMRT Link User Guide](#)):

1. Download the FASTA file containing the relevant barcode sequences from PacBio’s [Multiplexing](#) website or contact PacBio [Technical Support](#) to obtain a copy of the appropriate Barcode FASTA file. For example:
 - **HiFiViral_SARS-CoV-2_M13barcodes** FASTA file contains a list of 32 Forward and 32 Reverse M13 barcodes for use with the Barcoded M13 Primer Plate included in HiFiViral SARS-CoV-2 Kit ([102-132-000](#))

EXAMPLE FASTA FILE CONTAINING A LIST OF FORWARD AND REVERSE M13 BARCODES

```
>M13_bc1005_F
CACTCGACTCTCGCGTGAAAACGACGGCCAGT
>M13_bc1006_F
CATATATATCAGCTGTGAAAACGACGGCCAGT
>M13_bc1007_F
TCTGTATCTCTATGTGTTAAAACGACGGCCAGT
>M13_bc1008_F
ACAGTCGAGCGCTGCGGTA AAAACGACGGCCAGT
>M13_bc1009_F
ACACACGCGAGACAGAGTAAAACGACGGCCAGT
>M13_bc1010_F
ACGCGCTATCTCAGAGGTA AAAACGACGGCCAGT
>M13_bc1011_F
CTATACGTATATCTATGTA AAAACGACGGCCAGT
>M13_bc1012_F
```

```
1GTATCGCTCTATGTA AAAACGACGGCCAGT
>M13_bc1075_R
TAGAGAGCGTCGCGTGCAGGAAACAGCTATGAC
>M13_bc1076_R
GTGCACTCGCGCTCTCCAGGAAACAGCTATGAC
>M13_bc1077_R
TATCTCTCGAGTCGCGCAGGAAACAGCTATGAC
>M13_bc1078_R
CTCACACATACACGTCAGGAAACAGCTATGAC
>M13_bc1079_R
ATAGTACACTCTGTGTCCAGGAAACAGCTATGAC
>M13_bc1082_R
GTGACACACAGAGCACAGGAAACAGCTATGAC
```


IMPORTING THE BARCODE FASTA FILE INTO SMRT LINK FOR AUTOMATED DEMULTIPLEXING OF HiFiViral SARS-Cov-2 LIBRARY SAMPLES (CONT.)

2. Import the desired FASTA file into SMRT Link.

i. On the SMRT Link Home Page, select **Data Management**.

ii. Click **Import Data** and follow the steps below:

- 1 Specify whether to import data from the **SMRT Link Server**, or from a **Local File System**. (**Note:** Only references and barcodes are available if you select **Local File System**.)
- 2 Select the data type to import: **Barcodes** – FASTA (.fa or .fasta), XML (.barcodeset.xml), or ZIP files containing barcodes.
- 3 Navigate to the appropriate file and click **Import**. The selected barcode file is imported and becomes available for viewing in the SMRT Link Data Management module home screen.



The screenshot shows the 'Import Data' interface in the SMRT Link Data Management module. It is divided into three numbered steps:

- 1. Import from:** Two radio buttons are visible: 'SMRT Link Server' and 'Local File System'. The 'Local File System' option is selected.
- 2. Select File:** A dropdown menu for 'Data Type' is open, showing the following options: 'Barcodes (FASTA)' (highlighted in blue), 'Barcodes (ZIP)', 'References (FASTA)', and 'References (ZIP)'.
- 3. Browse Local Files:** A text input field contains 'Select File: (.fa, .fasta)'. To the right are 'Browse' and 'Import' buttons.

A 'Cancel' button is located at the bottom left of the interface.

SMRT LINK RUN DESIGN SETUP PROCEDURE FOR HiFiViral SARS-CoV-2 LIBRARY SAMPLES

A. Specifying Sample Information and Movie Collection Parameters

- 1 Under Application Type, select '**Viral Sequencing / HiFiViral SARS-CoV-2**'
 - Verify all default values in auto-filled sample information fields match the recommended values shown in **Quick Reference Card – Loading and Pre-Extension Time Recommendations for the Sequel II/IIe Systems** ([101-769-100](https://www.pacb.com/support/quick-reference-cards/sequel-ii-ii-e-quick-reference-card/)) for preparing HiFiViral samples for sequencing
- 2 Enter a Well Sample Name for your library sample
- 3 We recommend using a starting on-plate concentration (OPLC) = 200 pM and adjusting higher or lower if needed to achieve optimal *P1* loading

▼ SAMPLE 1: HiFiViral_SARS-CoV-2_Library_01, A01, 8 hour movie, 800 bp insert

Import from Sample Setup Select Sample

1 Application Required HiFiViral SARS-CoV-2

2 Well Sample Name Required HiFiViral_SARS-CoV-2_Library_01

Bio Sample Name Required [Empty]

Sample Comment [Empty]

Sample Well A01

Template Prep Kit Required SMRTbell® Express Template Prep Kit 2.0

Binding Kit Required Sequel® II Binding Kit 2.1

Sequencing Kit Required Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex Sequel® II DNA Internal Control 1.0

Insert Size (bp) Required 800

Recommended Concentration on Plate (pM) 100-300 pM

3 On-Plate Loading Concentration (pM) Required 200

Movie Time per SMRT Cell (hours) 8

Use Pre-Extension YES NO

Generate HiFi Reads ON INSTRUMENT IN SMRT LINK DO NOT GENERATE

SMRT LINK RUN DESIGN SETUP PROCEDURE FOR HiFiViral SARS-CoV-2 LIBRARY SAMPLES (CONT.)

B. Enabling Automated SARS-CoV-2 Data Analysis in SMRT Link and Specifying Sample Barcoding Information

1. To enable automated SARS-CoV-2 data analysis in SMRT Link:

- 1 Select YES for 'Automatic Launch of SARS-CoV-2 Analysis'
- 2 Enter an Analysis Name

Auto Analysis

1 Automatic Launch of SARS-CoV-2 Analysis YES NO

2 Analysis Name Required

2. Under Barcoded Sample Options, the following options are automatically specified if *HiFiViral SARS-CoV-2* is selected for Application Type:

- 3 Sample is Barcoded: **Yes**
- 4 Barcode Set: **HiFiViral_SARS-CoV-2_M13barcodes**
- 5 Same Barcodes on Both Ends of Sequence: **No**

Barcoded Sample Options

3 Sample Is Barcoded YES NO

4 Barcode Set Required

5 Same Barcodes on Both Ends of Sequence ? YES NO

SMRT LINK RUN DESIGN SETUP PROCEDURE FOR HiFiViral SARS-Cov-2 LIBRARY SAMPLES (CONT.)

3. Specify Barcode assignments and Bio Sample Names as follows:

- 6 Under **Assign Bio Sample Names to Barcodes**: Click **From a File**, then click **Download File**.
- 7 Edit the file and enter the **biological sample name**, **Plate ID** and **Plate Well** associated with each unique forward + reverse barcode pair listed in the first column; then save the file.
 - Delete entire rows of barcodes not used
 - Allowed characters*: Alphanumeric; dot; underscore; hyphen. Other characters will be automatically removed.
- 8 Browse for the Barcoded Sample File you just edited and click on Open.
- 9 You will see 'Upload was successful' appear on the line below, assuming the file is formatted correctly.

Assign Bio Sample Names to Barcodes Required Interactively From a File **6**

Autofilled Barcoded Sample File ? Download File

7

Barcode	Bio Sample Name	assayPlateID	assayPlateWell
M13_bc1002_F--M13_bc1050_R	HiFiViral_SARS-CoV-2_Sample_1	A	A01
M13_bc1002_F--M13_bc1051_R	HiFiViral_SARS-CoV-2_Sample_2	A	B01
M13_bc1002_F--M13_bc1052_R	HiFiViral_SARS-CoV-2_Sample_3	A	C01
M13_bc1002_F--M13_bc1053_R	HiFiViral_SARS-CoV-2_Sample_4	A	D01
M13_bc1017_F--M13_bc1070_R	HiFiViral_SARS-CoV-2_Sample_381	D	F11
M13_bc1017_F--M13_bc1071_R	HiFiViral_SARS-CoV-2_Sample_382	D	F12
M13_bc1017_F--M13_bc1072_R	HiFiViral_SARS-CoV-2_Sample_383	D	G12
M13_bc1017_F--M13_bc1073_R	HiFiViral_SARS-CoV-2_Sample_384	D	H12

8

Barcoded Sample Name File Required Browse

9

Barcoded Sample Name File Required Browse

Upload was successful

Refer to “Working with Barcoded Data” section in the [SMRT Link User Guide](#) for further details on how to specify barcode setup and sample name information in a Run Design



Multiplexed SARS-CoV-2 Data Analysis Recommendations

USE SMRT LINK TO EASILY ANALYZE MULTIPLEXED HiFi DATA FROM SARS-CoV-2 SURVEILLANCE SAMPLES

Analyze HiFiViral SARS-CoV-2 HiFi Data Using SMRT Link* by Creating an **Auto Analysis** in Run Design or by Performing a **Manual Analysis** in SMRT Analysis

Creating an Auto Analysis in Run Design

- **HiFiViral SARS-CoV-2 Analysis Application** can be run using the **Auto Analysis** feature available in SMRT Link Run Design
- This optional Run Design feature allows users to **automatically** complete all necessary analysis steps immediately after sequencing on the Sequel II and Ile Systems **without manual intervention**
- HiFiViral Auto Analysis workflow **automatically** launches CCS Analysis, Demultiplex Barcodes, and HiFiViral SARS-CoV-2 Analysis.


Performing a Manual Analysis in SMRT Analysis

- **HiFiViral SARS-CoV-2 Analysis Application** can also be run by performing a **manual analysis** in SMRT Link SMRT Analysis
- This process requires users to **manually** prepare input data for the HiFiViral SARS-CoV-2 Analysis Application
- HiFiViral manual analysis workflow requires **manually** specifying CCS Analysis ('Generate HiFi Reads') in Run Design, and then **manually** launching Demultiplex Barcodes and HiFiViral SARS-CoV-2 Analysis applications in SMRT Analysis

HiFiViral SARS-CoV-2 Auto Analysis Workflow



HiFiViral SARS-CoV-2 Manual Analysis Workflow



* Analysis is supported for samples isolated from individual humans and has not been designed or validated for use with other sample types (e.g., wastewater samples).

HiFiViral SARS-CoV-2 ANALYSIS SETUP – AUTO ANALYSIS

How to Use SMRT Link Run Design to Create an Auto Analysis

A. Specify Auto Analysis in Run Design

- 1 Under **Auto Analysis**, select YES for 'Automatic Launch of SARS-CoV-2 Analysis'
- 2 Enter an **Analysis Name**

Auto Analysis

1 Automatic Launch of SARS-CoV-2 Analysis YES NO

2 Analysis Name Required


B. Specify Barcoded Sample Options

Under **Barcoded Sample Options**, the following options are automatically specified if *HiFiViral SARS-CoV-2* is selected for Application Type in Run Design:

- 1 Sample is Barcoded: **Yes**
- 2 Barcode Set: **HiFiViral_SARS-CoV-2_M13barcodes**
- 3 Same Barcodes on Both Ends of Sequence: **No**

Barcoded Sample Options

1 Sample Is Barcoded YES NO

2 Barcode Set Required 

3 Same Barcodes on Both Ends of Sequence ? YES NO

How to Use SMRT Link Run Design to Create an Auto Analysis (Cont.)

- 4 Under **Assign Bio Sample Names to Barcodes**: Click **From a File**, then click **Download File**.
- 5 Edit the file and enter the **biological sample name**, **Plate ID** and **Plate Well** associated with each unique forward + reverse barcode pair listed in the first column; then save the file.
 - Delete entire rows of barcodes not used
 - Allowed characters*: Alphanumeric; dot; underscore; hyphen. Other characters will be automatically removed.
- 6 Browse for the Barcoded Sample File you just edited and click on Open.
- 7 You will see 'Upload was successful' appear on the line below, assuming the file is formatted correctly.

Assign Bio Sample Names to Barcodes Required Interactively **From a File** 4

Autofilled Barcoded Sample File ? **Download File**

5

Barcode	Bio Sample Name	assayPlateID	assayPlateWell
M13_bc1002_F--M13_bc1050_R	HiFiViral_SARS-CoV-2_Sample_1	A	A01
M13_bc1002_F--M13_bc1051_R	HiFiViral_SARS-CoV-2_Sample_2	A	B01
M13_bc1002_F--M13_bc1052_R	HiFiViral_SARS-CoV-2_Sample_3	A	C01
M13_bc1002_F--M13_bc1053_R	HiFiViral_SARS-CoV-2_Sample_4	A	D01
M13_bc1017_F--M13_bc1073_R	HiFiViral_SARS-CoV-2_Sample_381	D	F11
M13_bc1017_F--M13_bc1071_R	HiFiViral_SARS-CoV-2_Sample_382	D	F12
M13_bc1017_F--M13_bc1072_R	HiFiViral_SARS-CoV-2_Sample_383	D	G12
M13_bc1017_F--M13_bc1073_R	HiFiViral_SARS-CoV-2_Sample_384	D	H12

6

Barcoded Sample Name File Required Choose file **Browse**

7

Barcoded Sample Name File Required Barcode_Names_HiFiViral_384_Samples.csv **Browse**

Upload was successful


Refer to “**Working with Barcoded Data**” section in the **SMRT Link User Guide** for further details on how to specify barcode setup and sample name information in a Run Design

HiFiViral SARS-CoV-2 ANALYSIS SETUP – MANUAL ANALYSIS

How to Use SMRT Link SMRT Analysis to Perform a Manual Analysis

A. Prepare Input Data for the HiFiViral SARS-CoV-2 Analysis Application by Running Demultiplex Barcodes

- 1 In SMRT Analysis, select the SMRT Link **Demultiplex Barcodes** application, where the input to that application are HiFi Reads. (If HiFi Reads have not already been generated on the instrument, run CCS Analysis first.)
- 2 Barcode Set: Select **HiFiViral_SARS-CoV-2_M13barcodes**
- 3 Barcodes on Both Ends of Sequence: Select **No**



SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Demultiplex Barcodes

Import Analysis Settings Export

Associated Inputs

Barcode Set Required

HiFiViral_SARS-CoV-2_M13barcodes

Same Barcodes on Both Ends of Sequence

YES NO

How to Use SMRT Link SMRT Analysis to Perform a Manual Analysis (Cont.)

4 Under **Assign Bio Sample Names to Barcodes**: Click **From a File**, then click **Download File**.

5 Edit the file and enter the **biological sample name** associated with each unique forward + reverse barcode pair listed in the first column; then save the file.

- Delete entire rows of barcodes not used
- Allowed characters*: Alphanumeric; dot; underscore; hyphen. Other characters will be automatically removed.

6 Browse for the Barcoded Sample File you just edited and click on Open.

7 You will see 'Upload was successful' appear on the line below, assuming the file is formatted correctly.

8 Enter a Name for the Demultiplexed Output Data Set.

Assign Bio Sample Names to Barcodes Required ? Interactively From a File 4

Autofilled Barcoded Sample File ?

5

Barcode	Bio Sample Name
M13_bc1002_F--M13_bc1050_R	HiFiViral_SARS-CoV-2_Sample_1
M13_bc1002_F--M13_bc1051_R	HiFiViral_SARS-CoV-2_Sample_2
M13_bc1002_F--M13_bc1052_R	HiFiViral_SARS-CoV-2_Sample_3
M13_bc1002_F--M13_bc1053_R	HiFiViral_SARS-CoV-2_Sample_4
M13_bc1017_F--M13_bc1070_R	HiFiViral_SARS-CoV-2_Sample_381
M13_bc1017_F--M13_bc1071_R	HiFiViral_SARS-CoV-2_Sample_382
M13_bc1017_F--M13_bc1072_R	HiFiViral_SARS-CoV-2_Sample_383
M13_bc1017_F--M13_bc1073_R	HiFiViral_SARS-CoV-2_Sample_384

6

Barcoded Sample Name File ? Required

7

Barcoded Sample Name File ? Required Barcode_Names_HiFiViral_384_Samples.csv

8

Demultiplexed Output Data Set Name Required ?

HiFiViral_SARS-CoV-2_Sample_Plate_01_CCS_Demux

Refer to “**Working with Barcoded Data**” section in the [SMRT Link User Guide](#) for further details on how to specify barcode setup and sample name information in a Run Design

* **DO NOT include spaces** – Sample Names must be unique and will be truncated after any spaces.

How to Use SMRT Link SMRT Analysis to Perform a Manual Analysis (Cont.)

B. Set Up and Launch HiFiViral Analysis Application

1 After running the Demultiplex Barcodes application, create a new analysis using **SMRT Analysis > Create New Analysis**.

2 Name the analysis

3 Select **Data Types > HiFi Reads**.

4 Select all the demultiplex samples contained in the Data Set and choose **Analysis of Multiple Data Sets > One Analysis for All Data Sets**.

5 Under Analysis of Multiple Data Sets, specify '**One Analysis for All Data Sets**'

5 Click **Next**.

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Name Required: HiFiViral_SARS-CoV-2_Manual_Analysis_Demo

Analysis Type: AUTO ANALYSIS ANALYSIS

Data Type: HiFi Reads

Analysis of Multiple Data Sets: One Analysis for All Data Sets

Data Sets for selected Data Type displayed in table below.

Back


Members of HiFiViral_DataSet_96_Demux

	Data Set Details >				Sample Details			Run Data >	
<input checked="" type="checkbox"/>	Name	Well Sample Name	Run Name	Date Created	Created By	Bio Sample Name	Barcode Name	Total Length of Read	Instrument
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt117-83	M13_bc1014_F--M13_...	12,616,790	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt117-32	M13_bc1006_F--M13_...	56,575,682	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt117-09	M13_bc1002_F--M13_...	43,631,875	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt117-95	M13_bc1016_F--M13_...	141,207	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt114-29	M13_bc1006_F--M13_...	41,899,685	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt114-54	M13_bc1010_F--M13_...	27,074,587	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt114-27	M13_bc1006_F--M13_...	25,208,512	64011
<input checked="" type="checkbox"/>	HiFiViral_DataSet_D...	Twist 14-17-93well	20210917-Twist-Cr...	2021-09-20, 04:59:...	sizhang	Crt114-18	M13_bc1004_F--M13_...	22,770,872	64011

How to Use SMRT Link SMRT Analysis to Perform a Manual Analysis (Cont.)

B. Set Up and Launch HiFiViral Analysis Application (Cont.)

- 6 Select HiFiViral SARS-CoV-2 Analysis from the **Analysis Application** list.
- 7 Under **Associated Inputs**, SARS-CoV-2 Genome NC_045512.2 (the Wuhan reference genome) and Probe Sequences v1 are automatically loaded; advanced users may select a different reference or probe set if desired.
- 8 To generate the optional **Plate QC** graphical summary, click **Advanced Parameters** and load a CSV file using the provided template (assayPlateQC_template_4by96.csv) as a guide.



SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

HiFiViral SARS-CoV-2 Analysis

Import Analysis Settings Export

Associated Inputs

Reference Genome Required

SARS-CoV-2 Genome NC_045512.2

Probe sequences Required

HiFiViral SARS-CoV-2 Enrichment Probes

Advanced Parameters

How to Use SMRT Link SMRT Analysis to Perform a Manual Analysis (Cont.)

B. Set Up and Launch HiFiViral Analysis Application (Cont.)

9 Under **Advanced Parameters**, download the provided CSV template (`assayPlateQC_template_4by96.csv`) as a guide and edit the file.

Enter the **biological sample name**, **Plate ID** and **Plate Well** associated with each unique forward + reverse barcode pair listed in the first column; then save the file.

- Delete entire rows of barcodes not used
- Allowed characters*: Alphanumeric; dot; underscore; hyphen. Other characters will be automatically removed.

Browse for the Plate QC File you just edited and click on Open.

You will see 'Upload was successful' appear on the line below, assuming the file is formatted correctly.

10 Click **Start** to start the analysis.

Advanced Analysis Parameters

Plate QC CSV Minimum Base Coverage Minimum Variant Frequency

Minimum Barcode Score Advanced Processing Options Compute Settings

Barcode	Bio Sample Name	assayPlateID	assayPlateWell
M13_bc1002_F--M13_bc1050_R	HiFiViral_SARS-CoV-2_Sample_1	A	A01
M13_bc1002_F--M13_bc1051_R	HiFiViral_SARS-CoV-2_Sample_2	A	B01
M13_bc1002_F--M13_bc1052_R	HiFiViral_SARS-CoV-2_Sample_3	A	C01
M13_bc1002_F--M13_bc1053_R	HiFiViral_SARS-CoV-2_Sample_4	A	D01
M13_bc1002_F--M13_bc1054_R		A	E01
M13_bc1002_F--M13_bc1055_R		A	F01
M13_bc1002_F--M13_bc1056_R		A	G01
M13_bc1002_F--M13_bc1057_R		A	H01
M13_bc1002_F--M13_bc1058_R		A	A02
M13_bc1002_F--M13_bc1059_R		A	B02
M13_bc1002_F--M13_bc1060_R		A	C02
M13_bc1002_F--M13_bc1061_R	HiFiViral_SARS-CoV-2_Sample_12	A	D02
M13_bc1002_F--M13_bc1062_R	HiFiViral_SARS-CoV-2_Sample_13	A	E02

Plate QC CSV

Upload was successful

How to Use SMRT Link SMRT Analysis to Perform a Manual Analysis (Cont.)

Comparison of CSV Templates for Demultiplex Barcodes Analysis and HiFiViral SARS-CoV-2 Assay Plate QC Analysis

Demultiplex Barcodes

Barcode	Bio Sample Name
M13_bc1002_F--M13_bc1050_R	HiFiViral_SARS-CoV-2_Sample_1
M13_bc1002_F--M13_bc1051_R	HiFiViral_SARS-CoV-2_Sample_2
M13_bc1002_F--M13_bc1052_R	HiFiViral_SARS-CoV-2_Sample_3
M13_bc1002_F--M13_bc1053_R	HiFiViral_SARS-CoV-2_Sample_4
M13_bc1017_F--M13_bc1070_R	HiFiViral_SARS-CoV-2_Sample_381
M13_bc1017_F--M13_bc1071_R	HiFiViral_SARS-CoV-2_Sample_382
M13_bc1017_F--M13_bc1072_R	HiFiViral_SARS-CoV-2_Sample_383
M13_bc1017_F--M13_bc1073_R	HiFiViral_SARS-CoV-2_Sample_384

CSV Template contains two columns

HiFiViral SARS-CoV-2 Assay Plate QC

Barcode	Bio Sample Name	assayPlateID	assayPlateWell
M13_bc1002_F--M13_bc1050_R	HiFiViral_SARS-CoV-2_Sample_1	A	A01
M13_bc1002_F--M13_bc1051_R	HiFiViral_SARS-CoV-2_Sample_2	A	B01
M13_bc1002_F--M13_bc1052_R	HiFiViral_SARS-CoV-2_Sample_3	A	C01
M13_bc1002_F--M13_bc1053_R	HiFiViral_SARS-CoV-2_Sample_4	A	D01
M13_bc1017_F--M13_bc1070_R	HiFiViral_SARS-CoV-2_Sample_381	D	F01
M13_bc1017_F--M13_bc1071_R	HiFiViral_SARS-CoV-2_Sample_382	D	F12
M13_bc1017_F--M13_bc1072_R	HiFiViral_SARS-CoV-2_Sample_383	D	G12
M13_bc1017_F--M13_bc1073_R	HiFiViral_SARS-CoV-2_Sample_384	D	H12

CSV Template contains four columns

When editing CSV templates for Demultiplex Barcodes analysis and HiFiViral SARS-CoV-2 Assay Plate QC analysis:

- ❑ Delete entire rows of barcodes not used
- ❑ Allowed characters*: Alphanumeric; dot; underscore; hyphen. Other characters will be automatically removed.
 - **DO NOT** include spaces – Sample Names must be unique and will be truncated after any spaces.

HiFiViral SARS-CoV-2 ANALYSIS WORKFLOW

SMRT Link HiFiViral SARS-CoV-2 Auto Analysis* Workflow Algorithm Descriptions



- 1. Demultiplex barcodes** using the `lima` tool, where the input to that application are HiFi Reads HiFi (\geq Q20 CCS) Reads (BAM format).
- 2. Process the reads to trim the probe arm sequences** using the `mimux` tool.
- 3. Align the reads** to the reference genome using `pbmm2`.
- 4. Call and filter variants** using `bcftools`, generating the raw variant calls in VCF file format. Filtering in this step removes low-quality calls (less than Q20), and normalizes indels.
- 5. Filter low-frequency variants** using `vcfcons` and generate a consensus sequence by injecting variants into the reference genome. At each position, a variant is called only if both the base coverage exceeds the minimum base coverage threshold and the fraction of reads that support this variant is above the minimum variant frequency threshold.



* The SMRT Link Demultiplex Barcodes and HiFiViral SARS-CoV-2 Analysis Applications must each be launched manually if Auto Analysis is not specified in Run Design when setting up a sequencing run on Sequel II or Ite Systems with HiFiViral SARS-CoV-2 Kit library samples. 49

HiFiViral SARS-CoV-2 ANALYSIS OUTPUTS

SMRT Link HiFiViral SARS-CoV-2 Analysis Application Outputs

- Per-sample analysis outputs include:
 - ❑ Consensus sequence (FASTA)
 - ❑ Variant calls (VCF)
 - ❑ HiFi Reads aligned to the reference (BAM)
 - ❑ Sample Summary table including: Count of variable sites, genome coverage, read coverage, and probability of multiple strains, and other metrics
 - ❑ Plot of HiFi Read coverage across the SARS-CoV-2 genome

Sample Summary


Bio Sample Name	Substitutions	Insertions	Deletions	Reads	Read Coverage	On-Target Rate	Multiple Strains (Probability)	Ns	Genome Coverage
Sample 1	36	0	3	12,964	288	99.99%	No (0.00)	156	99.47%
Sample 2	38	0	3	1,075	24	99.81%	No (0.00)	761	97.45%
Sample 3	40	0	3	2,289	51	99.91%	No (0.00)	219	99.26%

File Downloads

Edit Output File Name Prefix Example: analysis-Twist Bioscience Control 17-136917

File	Size	Type
All Samples, Probe Counts TSV	935 KB	zip
Sample Summary Table CSV	9 KB	csv
All Samples, Raw Variant Call VCF	267 KB	zip
All Samples, Consensus Sequence Aligned BAM	819 KB	zip
All Samples, HiFi Reads Mapped BAM	515 MB	zip
All Samples, Variant Call VCF	250 KB	zip
Analysis Log	737 KB	log
All Samples, Genome Coverage Plots	30 MB	zip



Genome Coverage



SMRT Link HiFiViral SARS-CoV-2 Analysis Application Outputs (Cont.)

- HiFiViral SARS-CoV-2 analysis application also outputs a graphical summary of performance across all samples in assay plate layout for Sample Plate QC evaluation

Lower HiFi read counts due to **evaporation-induced edge-effects** during viral enrichment



DOWNLOADING HiFiViral SARS-CoV-2 ANALYSIS RESULTS IN SMRT LINK V10.2

To download the HiFiViral SARS-CoV-2 analysis results, click on the File Downloads tab to download the desired output files.

- ▶ Analysis Overview
- ▶ Summary Report
- ▼ Data
- File Downloads**
- SMRT Link Log

File Downloads

Edit Output File Name Prefix **Example:**analysis-Twist_RNA_23_Ct29p8_rep2-45495

File	Size	Type
📁 All Samples, Probe Counts TSV	994 KB	zip
📁 Sample Summary Table CSV	9 KB	csv
📁 All Samples, Raw Variant Call VCF	244 KB	zip
📁 All Samples, Consensus Sequence Aligned BAM	791 KB	zip
📁 All Samples, HiFi Reads Mapped BAM	707 MB	zip
📁 All Samples, Variant Call VCF	206 KB	zip
📁 All Samples, Genome Coverage Plots	33 MB	zip
📁 All Samples, Consensus Sequence FASTA	701 KB	zip
📁 All Samples, HiFi Reads FASTQ	871 MB	zip
📁 Analysis Log	761 KB	log
📁 Analysis Log	25 KB	log

DOWNLOADING HiFiViral SARS-CoV-2 ANALYSIS RESULTS IN SMRT LINK V10.2 (CONT.)

analysis-Twist_RNA_23_Ct29p8_rep2-45495-samples.consensus.fasta	--	Folder
Twist_RNA_13_Ct19p1_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct19p1_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct19p1_rep3.consensus.fasta		
Twist_RNA_13_Ct21p9_rep1.consensus.fasta		
Twist_RNA_13_Ct21p9_rep2.consensus.fasta		
Twist_RNA_13_Ct21p9_rep3.consensus.fasta		
Twist_RNA_13_Ct22p6_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct22p6_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct22p6_rep3.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct24p4_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct24p4_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct24p4_rep3.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct26p2_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct26p2_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct26p2_rep3.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct28_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct28_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct28_rep3.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct29p8_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct29p8_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct29p8_rep3.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct31p5_rep1.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct31p5_rep2.consensus.fasta	30 KB	FASTA File
Twist_RNA_13_Ct31p5_rep3.consensus.fasta	30 KB	FASTA File
Twist_RNA_14_Ct19p1_rep1.consensus.fasta	30 KB	FASTA File

For each sample, HiFiViral analysis application outputs a single SARS-CoV-2 consensus sequence




Multiplexed SARS-CoV-2 Library Example Performance Data


EXAMPLE SEQUENCING PERFORMANCE FOR TWIST SYNTHETIC SARS-CoV-2 RNA CONTROLS [6 X 5 KB FRAGMENTS]

SMRTbell Library QC and Primary Sequencing Metrics for 96-Plex and 384-Plex Twist Control Samples


Library QC



Raw Data Report



CCS Analysis Report



	96-Plex	384-Plex
Yield of Pooled Barcoded PCR Products	2049 ng	12,400 ng
Pooled DNA Input for Library Construction	1000 ng	1000 ng
Final Yield of ProNex Bead Purified Library (%)	142 ng (14.2%)	408 ng (40.8%)

Twist 14 and 17 controls (were enriched with the HiFiViral SARS-CoV-2 Kit. Pooled barcoded PCR products were purified with 1.3X ProNex Beads and constructed into SMRTbell libraries with SMRTbell Express TPK 2.0.

	96-Plex	384-Plex
Raw Base Yield	145.6 Gb	139.4 Gb
Mean Polymerase Read Length	26.3 kb	25.1 kb
P0	18.9%	19.7%
P1	69.2%	69.4%
P2	11.9%	10.9%

200 pM on-plate concentration / 8-h movie time / No Pre-Extension Time / No Adaptive Loading

	96-Plex	384-Plex
HiFi Reads	3.6 M	3.5 M
HiFi Base Yield	2.8 Gb	2.8 Gb
Mean HiFi Read Length	772	788
Median HiFi Read Quality	QV60	QV60
HiFi Read Mean # of Passes	21	21

EXAMPLE SEQUENCING PERFORMANCE FOR TWIST SYNTHETIC SARS-CoV-2 RNA CONTROLS [6 X 5 KB FRAGMENTS] (CONT.)

HiFiViral SARS-CoV-2 Auto Analysis Outputs for 96-Plex Twist Control Samples

Summary Report


Value	Analysis Metric
96	Samples
93	Samples With Genome Coverage > 90%
92	Samples With Genome Coverage > 95%
0	Samples Failing Workflow

Sample Summary

Bio Sample Name	Substitutions	Insertions	Deletions	Reads	Read Coverage	On-Target Rate	Multiple Strains (Probability)	Ns	Genome Coverage
Ctrl17-96	0	0	0	5	0	100.00%	No (0.00)	29,903	0.00%
Ctrl17-31	32	1	1	55,235	1,197	99.99%	No (0.00)	616	97.94%
Ctrl_14-27	31	0	4	35,341	762	100.00%	No (0.00)	682	97.72%
Ctrl_14-03	30	0	4	9,362	177	100.00%	No (0.02)	1,556	94.79%

- 93 Positive Control samples showed ≥90% genome coverage (Blue and Green wells in Plate QC image)
- 3 Negative Control samples showed <90% genome coverage as expected (Yellow wells)

Plate QC



HiFiViral SARS-CoV-2 KIT DELIVERS ROBUST GENOME COVERAGE PERFORMANCE ACROSS VARIABLE INPUT QUANTITIES AND MULTIPLEX LEVELS

Example SARS-CoV-2 Genome Coverage Results Obtained for Twist Control Samples

Experimental Design

96-plex prepared with 4 Synthetic Twist RNA Controls at 8 input quantities in replicates of 3.

TWIST CONTROL	VARIANT	PART NUMBER
14	Alpha (B.1.1.7)	103907
15	Alpha (B.1.1.7)	103909
16	Beta (B.1.351)	104043
17	Gamma (P.1)	104044


RNA Input Quantity*

SAMPLE CT	COPY NUMBER
19	6 M
20	3 M
21	1 M
24	100,000
27	10,000
30	1,000
33	100
35	3

Input Quantity Input of RNA controls ranged from 6 million copies down to 3. Copy number is converted into Ct scale after Han *et al.* 2021.*

* Han M.S., et al. (2021). RT-PCR for SARS-CoV-2: quantitative versus qualitative. *The Lancet Infectious Disease* 21(2) p165.

96-Plex of Twist Control Samples



Samples with Ct ≤ 30 achieved complete genome coverage** with 1000 HiFi reads


- 4-fold HiFi Read depth required to output a consensus base
- $\sim 1,000$ mapped HiFi reads reliably yields $\geq 90\%$ genome coverage

** Complete = $\geq 90\%$ genome coverage


HiFiViral SARS-CoV-2 KIT DELIVERS ROBUST GENOME COVERAGE PERFORMANCE ACROSS VARIABLE INPUT QUANTITIES AND MULTIPLEX LEVELS (CONT.)

Example SARS-CoV-2 Genome Coverage Results Obtained for Surveillance Samples

96-plex of "Real" Samples for Surveillance



384-plex of Controls and Nasopharyngeal Extracts



Genome Completeness in Surveillance Samples

SAMPLE INPUT	NO. OF SAMPLES	> 90% GENOME COVERAGE
Known Ct	84	83%
Unknown Ct	9	44%
Twist Controls	2	100%
Negative Control	1	0


Genome Completeness in 384-plex

SAMPLE INPUT	NO. OF SAMPLES	> 90% GENOME COVERAGE
Controls (Ct<30)	216	90%
NP Extracts	144	85%

HiFiViral SARS-CoV-2 KIT ENABLES COMPREHENSIVE CHARACTERIZATION OF VARIANTS FOR SURVEILLANCE AND COVID-19 RESEARCH

SARS-CoV-2 Variant Calling Achieves High Precision and Recall for Characterization of SNVs and SVs



Recovery of Large Deletions in ORF7a



Deletions (87 – 271 bp) are called in VCF and consensus sequence.

SNV Calling & Strain Assignment for Controls in a 384-plex

CONTROL SAMPLE	NEXTCLADE ASSIGNMENT	COMPLETE GENOMES	PRECISION	RECALL	NEXTSTRAIN ACCURACY
Twist 01	19A	29	1	94.8%	100%
Twist 13	20C	24	1	99.7%	100%
Twist 14	20I (Alpha, V1)	25	1	99.9%	100%
Twist 15	20I (Alpha, V1)	24	1	99.9%	100%
Twist 16	20H (Beta, V2)	24	1	100%	100%
Twist 17	20J (Gamma, V3)	24	1	100%	100%
Twist 23	21A (Delta)	24	99.1%	99.4%	100%



Example visualizations of HiFi reads spanning around large deletions.


- HiFi reads can detect SNVs and SVs with high precision and recall for accurate SARS-CoV-2 strain assignment

HiFiViral SARS-CoV-2 KIT ENABLES DETECTION OF MINOR VARIANTS AND MULTIPLE STRAINS* IN THE SAME SAMPLE


Mixed Control Experiment

- Titrated mixed controls
- Minor frequency: 1% to 50%
- Binomial model for multi-strain detection*
- Achieve $P > 95\%$ at $>20\%$ minor frequency**

Multi-Strain Calling Performance for Mixed Controls



Detection of Minor Variants in Surveillance Samples



Possible Sources of Multiple Strains in Sample

- Sample contamination, lab error, infection with multiple strains
- We recommend users confirm presence of multiple strains with additional experiments

* Multi-strain detection is supported for samples with Ct < 26

** Power of detection increases with more variable sites.



Technical Documentation & Applications Support Resources

TECHNICAL RESOURCES FOR SARS-CoV-2 LIBRARY PREPARATION, SEQUENCING & DATA ANALYSIS

Visit PacBio's [COVID-19 Sequencing Tools and Resources Website](#) for HiFiViral SARS-CoV-2 Workflow Updates and Other Resources

Sample Preparation Literature

- Procedure & Checklist – PacBio HiFiViral High-Throughput Multiplexing for Full-Viral Genome Sequencing of SARS-CoV-2 ([102-188-800](#))
- Quick Reference Card – Loading and Pre-extension Recommendations for the Sequel II/IIe Systems ([101-769-100](#))
- Overview – Sequel Systems Application Options and Sequencing Recommendations ([101-851-300](#))
- Application Brief: HiFiViral Full-Viral Genome Sequencing – Best practices ([BP110-111121](#))
- Application Note: HiFiViral Full-Viral Genome Sequencing (102-194-700) [*Coming Soon*]
- Technical Overview: Multiplexed Library Preparation for Full-Viral Genome Sequencing Using HiFiViral SARS-CoV-2 Kit ([102-205-300](#))

FAQ

- HiFiViral SARS-CoV-2 Kit FAQ [[Link](#)]

Posters, Videos & Webinars

- PacBio HiFiViral SARS-CoV-2 Kit Product Overview Video (2021) [[Link](#)]
- SFAF Poster (2021): HiFiViral SARS-CoV-2: A kitted solution for genome surveillance that is robust across sample input quantities and new variants [[Link](#)]
- ASHG Webinar (2021): HiFiViral SARS-CoV-2 Kit: A differentiate solution for surveillance by sequencing [[Link](#)]



PacBio HiFiViral SARS-CoV-2 Kit Product Overview Video (2021) [[Link](#)]

TECHNICAL RESOURCES FOR SARS-CoV-2 LIBRARY PREPARATION, SEQUENCING & DATA ANALYSIS (CONT.)

Ordering Information

CONSUMABLE PRODUCT	PART NUMBER
HiFiViral SARS-CoV-2 Kit (384 rxn)	102-132-000
SMRTbell Express Template Prep Kit 2.0 (18 rxn)	100-938-900
SMRT Cell 8M Tray	101-389-001
Sequel II Binding Kit 2.1 and Internal Control 1.0 (24 rxn)	101-843-000
Sequel II Sequencing Kit 2.0 (4 rxn)	101-820-200
SMRTbell Enzyme Cleanup Kit 2.0 (10 rxn)	101-932-600



APPENDIX 1: RNA Isolation Kit Options for Full-Viral Genome Sequencing of SARS-CoV-2

RNA SAMPLE EXTRACTION KIT OPTIONS FOR FULL-VIRAL GENOME SEQUENCING OF SARS-CoV-2

Note: The products below have not been tested or validated by PacBio but are listed here as examples of third-party kits used by other PacBio customers for isolating SARS-CoV-2 RNA samples for multiplexed SMRTbell library preparation

VENDOR	RNA ISOLATION KIT PRODUCT	AUTOMATION PLATFORM
Thermo Fisher Scientific	MagMAX Viral and Pathogen Nucleic Acid Isolation Kit [Link]	KingFisher Flex System
Roche Molecular Systems	MagNA Pure 96 DNA and Viral NA Small Volume Kit [Link]	Roche MagNA Pure-96 (MP6)




APPENDIX 2: Guidance on Workflow Automation For Multiplexed Library SARS-CoV-2 Library Preparation

WORKFLOW AUTOMATION OPTIONS FOR HIGH-THROUGHPUT MULTIPLEXED HiFiViral SARS-CoV-2 SAMPLE PREPARATION

Interested in automating your HiFiViral SARS-CoV-2 sample preparation workflow to achieve higher throughput? Please [contact PacBio Support](#) or your local Field Applications Scientist to discuss your needs.



Agilent Bravo Liquid Handler



Sciclone G3 NGSx Workstation



Biomek 4000 Workstation



Hamilton Microlab VANTAGE Liquid Handler



Tecan Infinite F-Series Plate Reader




Custom Liquid Handler

Key Considerations for Workflow Automation


- Liquid handler capabilities, including:
 - Small volume ($\geq 2 \mu\text{L}$) and large volume ($\geq 200 \mu\text{L}$) transfers
 - Magnetic plate blocks for bead-based purification and buffer exchanges
 - Integrated heating / cooling temperature control
- Microplate reader for high-throughput DNA concentration QC

RECOMMENDED STEPS TO AUTOMATE FOR VIRAL GENOME ENRICHMENT WORKFLOW USING HiFiViral SARS-CoV-2 KIT



* Off-deck processing steps shown may be moved on-deck depending on liquid handler capabilities

RECOMMENDED STEPS TO AUTOMATE FOR HiFiViral SARS-CoV-2 SMRTBELL LIBRARY CONSTRUCTION WORKFLOW





PACBIO®

www.pacb.com

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2021 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.